

#### ZIRCONIUM MILLING TECHNOLOGY

Move the world with your hands



### The Inventor.

Dear colleagues,

Please allow me to look back on the development history of our milling system.

Five years ago, I had the idea of elaborating zirconium by the means of a manual milling unit. A lot of patience and inventiveness was needed from the first idea to the actual realization of the milling unit ready for sale. I started with a wooden model and tried out several possibilities, before I developed a working prototype. After two years and 15,000 milled units, the first unit was ready for sale. Thanks to international teamwork we were able to sell our system in 60 different countries.



This brochure will tell you more about how the machine works. It would be a pleasure for me, if after having read it, you feel like trying out my system.

I am looking forward to your reaction!

Sincerely,

finico Stejen

MANPOWER	Move the world with your hands
	This is the motto of our milling system. It provides the technician with the possibility to make high-quality units with his own hands.



**Zirconium** is one of the oldest and most abundant elements in the terrestrial crust and it is the basis for Zirconium oxid (yttrium stabilized Zirconium dioxid). This first-class-performance material was successfully used for artificial limbs and joints in the medical field in the last decades and now it is also available to the dental industry. Due to its excellent biological characteristics, Zirconium is nowadays the preferred material for dental restorations. Its use in the dental field is on the increase since the 1990s. It is assumed that 15,000 to 20,000 units are made in zirconium every day. Laboratory tests revealed that zirconium bridges and metal ceramic bridges on the basis of precious metal alloy show the same fracture resistance.

«ZIRCONIAthe ceramic material of the future»

#### MATERIAL AND WORK EXPENDITURE

Milling bur consumption:	approx. 80-100 units per bur		
Milling time per unit:	approx. 10-15 min. per unit		
Time for modelling:	approx. 10-15 min. per unit		
Sinterization cycle:	approx. 8 hours		
Maintenance:	Cleaning and lubrication		
Wear and tear:	Wearing parts can be easily exchanged		
Consumption furnace:	900 Watt		
Total costs per unit, incl. all costs also ceramic:	15€		
Preparation:	All kinds of preparation possible		

IMPROVE YOUR MILLING TECHNIQUE	Our product range includes special instruments for milling all kinds of abutments and other demanding works.	3
		And a second



Designer Zirkograph

## Zirconium milling technology.

Zirkonzahn Ltd offers individual solutions of manual milling technology.

Technical solutions - for simple as well as highly complex and demanding constructions.

A wide range of accessory material rounds off the system.

#### ICE ZIRCONIA TRANSLUCENT

COMPOSITION			Specification
$\operatorname{Zr}O_2(+HfO_2)$	%	:	Main component
Y <sub>2</sub> 0 <sub>3</sub>	%	:	4.95 ~ 5.26
Al <sub>2</sub> 0 <sub>3</sub>	%	:	0.15 ~ 0.35
Si0 <sub>2</sub>	%	:	Max. 0.02
$\operatorname{Fe}_2 0_3$	%	:	Max. 0.01
Na <sub>2</sub> 0	%	:	Max. 0.04
Density (g/cm <sup>3</sup> ) sintered			6,05
Hardness (HV10)			>1250
Weibull modulus			>10
Flexural strength R.T.(MPa)		>1200 (MPa)	

#### ICE ZIRCONIA TRANSLUCENT AND ZIRCONIA PRETTAU

Both ICE Zirconia materials can be used for crowns and bridges. Due to the high level of translucency the Zirconia Prettau is especially suitable for full zirconia bridges.

#### SOME COMPONENTS OF THE SYSTEM

- 1. Milling unit "Designer Zirkograph"
- 2. Sinter furnace "Zirkonofen 600"
- 3. Infrared Predying lamp
- **4.** 5th axis
- 5. Assortment of ceramics
- **6.** Assortment of stains
- 7. Zirconia blanks
- 8. Colour liquid (for the colorization of zirconia)



















## Mock-up frame construction.



No limits on the type of preperations.

4.



Round off sharp corners – block out undercuts (Vaseline)



Apply Rigid T



Apply Rigid T stopping 1 mm above the margin and cure





Glue



Finalize crowns - do not elaborate



Insert bridge pontic (light cure tray material)



Cut bridge pontic

Reconnect bridge pontic (without plaster basis)



Resin template (Frame)

Mark perimeter on template Marked perimeters



Draw connectors

Cut with tungsten bur



Cut template

Fix bridge into template



Double check fixed bridge on model

## Milling procedure.



Start the milling process with the 4L bur. Use smaller instruments for fine and final milling.



Glue Zirconia block into position at both ends

Mill the outer form by using the bur 4L



Do not push too much during milling

Remove surplus material of the outer margin



Mill down to margin level

Internal milling of crowns







Precise internal milling and smoothing with 1L bur



Milled unit - ready for disconnection



Disconnected bridge with remaining support base – ready for immersion into Colour Liquid



Dip into Colour Liquid for 5 seconds





Put the objects under the infrared predying lamp. Overnight sinterization – firing cycle approx. 8 hours.

## Stratification.



"Much dentin – little enamel" Zilio Aldo, Venice



Elaborate zirconium frame structure with zirconium dioxid stones



Blast structure with aluminium oxyd at bar



Apply washbrand with dentin opaquer (high fluorescence)



Increase temperature 100°, holding time 2 min.



Surface should be as smooth as glass



Apply dentin with 50% of dentin orange



Apply dentin



Indentation of dentin in incisal areas



Apply thick coat of transpa 3 at the sides



Apply transpa blue into the approximal space



Apply the corresponding enamel mass Apply transpa masses in forms of stripes





First intermediate firing



Apply only dentin on the cervical area Add enamel





Apply T3 on incisal and approximal area



Bridge - after final firing



Zilio Aldo, Venice

### Cases.









### Cases.











South Tyrol – the magic of diversity

Nice sociable people and good food make your stay unforgettable





Enjoy with us South Tyrol's nature and calmness



### Introductory courses.

Every week we offer several introductory courses on the zirconia milling systems in our education centre in Bruneck (South Tyrol).



#### PROGRAM:

- Introduction to the milling concept
- Manufacture of resin frame work
- Tension-free frame assembly
- Milling zirconia
- Frame refinement at presintered stage
- Colouring
- Sinter fire over night
- Fitting and frame preparation for bonding
- Other hints and tricks

#### COURSE DURATION:

Day 1: 9 am to 6 pm (approx.) Day 2: 9 am to 12 pm (approx.)

Cost per person: 190,00 € (excl. tax) Maximum participants: 12

## REGISTRATION AND FURTHER INFORMATION:

www.zirkonzahn.com TEL: +39 0474 066 670





### ZIRCONIUM MILLING TECHNOLOGY

**Zirkonzahn World Wide** - An der Ahr 7 - 39030 Gais/South Tirol (Italy) T +39 0474 066 680 - F +39 0474 066 661 - www.zirkonzahn.com - info@zirkonzahn.com

