Gebrauchsanweisung - Informazione d'uso Instruction manual - Instrucciones de uso Mode d'emploi

Deutsch

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Date of update: 01.03.2012

Intended use:

After being modelled and adapted, the product is inserted into the human body by specialised technicians. It is used to model dental prosthesis substituting worn out or defective anatomic parts. This product is not destined to manufacture implants.

Manual for: ICE Zirkon and Prettau Zirkon

1. Fix the zirconia block in the specific clamping device.
2. Effectuate the milling process.
3. Remove milled frame work from block and pre-trim with diamond disc
4. Finish with tungsten cross-cut burs
5. De-dust frame with air gun
6. Colour frame with COLOUR LIQUID (see manual for COLOUR LIQUID)
7. Final-dry frame under infrared lamp for 45 minutes
8. Load furnace
9. Shut gate and select program. (Observe chart of sintering parameters)
10. Remove sintered framework from furnace at room temperature.
11. Check frame on working model for seat and margin fit. Trim with water cooled air turbine if required.
12. Blast using 120 micron Alu-oxide at pressure 4 bar.
13. Clean frame with brush using water. Do not steam-clean.
14. Frame is now ready for veneering

Chart of sintering parameters - final temperature

Material	Final temp.
Prettau® Zirconia	1600°C
Prettau® Zirconia Pre-Coloured	1450°C
Other types of Zirconia (e. g. ICE Zirkon Transluzent)	1500°C

Chart of sintering parameters - heating rates (valid for all materials)

Weight (incl. sintering foot/dental unit)	Heating rate	Holding time
0-2 g	max. 14°C/min.	1,5 h
2-4 g	max. 8°C/min.	2 h
over 4 g	max. 5°C/min.	2 h

Advices:

The product must be used only by specifically trained personnel. Should intolerance to the product be suspected, use it only after an allergy screening test has been performed with issue of medical statement of non-allergy conditions. This instruction manual is state of the art and is provided based on our current levels of knowledge. The user will be held fully responsible for any operations entailing the use of the product. The manufacturer is not liable for any erroneous results, since he has no influence on any processing procedure. Any claims for damages therefore shall concern exclusively the commercial value of our products.

The dental technician has to make sure that the dentist has received all of the necessary information for the correct use of the product.

Safety instructions:

Ensure adequate ventilation.

For the porcelain veneering of bridges of various elements you have to observe the following parameters during the heating- and the cooling-phase:

- 0 - 1 g per dental-unit 55° / min.
- 1 - 2 g per dental-unit 45° / min.
- 2 - 3 g per dental unit 35° / min.
- more than 3 g per dental unit 25° / min.

Wear protective glasses.

Protect respiratory system:

Wear suitable face mask to prevent „fine dust“ inhalation. Avoid „fine dust“ air pollution by using the air blaster in short bursts only rather than a continuous flow.

Storage and Handling:

Keep container dry. Keep at temperatures: < 28°C / < 82°F
CAUTION: US Federal law restricts this device to sale by or on the order of a licensed dentist.

Information

Producer: Zirkonzahn S.r.l.
An der Ahr 7 - 39030 Gais - Italy
www.zirkonzahn.com
Supplier: Zirkonzahn S.r.l.
An der Ahr 7 - 39030 Gais - Italy
www.zirkonzahn.com
Emergency Phone: +39 0474 066 660
(Mon. - Fri. 8 am - 6 pm)

ICE ZIRKON and PRETTAU ZIRKON have especially been designed for the use with ICE ZIRKON MALFARBEN, ICE ZIRKON KERAMIK and COLOUR LIQUID.

N.B.: All medical products are non-returnable.

Gebrauchsanweisung - Informazione d'uso
Instruction manual - Instrucciones de uso
Mode d'emploi

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