

FL-BOND II

Self-Etching Fluoride
Releasing Bonding System



PRG
Technology
Giomer



FL-BOND II

FL-BOND II is a self-etching fluoride releasing two-step adhesive system which features an excellent bond to both enamel and dentin and provides a secure marginal seal.

The system is characterized by its excellent handling and enhanced biocompatible properties.

Step 1 - Priming

The FL-BOND II Primer contains a new effective adhesive promoting monomer and new photo initiator and is free of HEMA and acetone to considerably minimize odour, post operative sensitivity and gingival white blanching.

The FL-BOND II self-etching Primer prepares the enamel and dentin surfaces while eliminating the need for separate etching with phosphoric acid.

Step 2 - Bonding

The FL-BOND II Bonding Agent incorporates fluoride containing S-PRG filler (Surface Pre-Reacted Glass-ionomer) which continuously provide fluoride to the tooth structure.

Because of its high filler loading FL-BOND II provides exceptional resistance against abrasion and effectively avoids a marginal washout along with durable sealing.




Indications


- Adhesive bonding of any composite to tooth structure
- Pre-treatment of the tooth structure prior to the application of resin cements
- Repair of restorations

Comfortable and quick handling


All it takes is two simple steps and about 35 seconds, because technique-sensitive or laborious steps such as etching, shaking, mixing or moisture surveillance are no longer necessary.




1 Dispense Primer, apply on enamel and dentin and leave undisturbed for 10 seconds



2 Dry with oil-free air (Do not rinse with water)




3 Dispense Bonding Agent and apply an even layer on the entire cavity



4 Light-cure: 10 seconds with halogen light or 5 seconds with LED light


Optimal viscosity and adhesion

The moderately viscous consistency of FL-BOND II allows a convenient application and is particularly stable with sloping cavity walls. This guarantees an even film thickness at all preparation surfaces and a stable bond strength.




Radiopacity for accurate diagnosis

For a reliable diagnosis FL-BOND II offers a remarkable radiopacity which provides an exceptional visualization of the restoration and the bonding layer on the X-ray.



Durable protection


FL-BOND II Bonding Agent contains 40% fluoride releasing and recharging S-PRG filler. This offers a permanent additional protection against secondary caries due to an optimized remineralization of the tooth structure.¹



Fluoride release and recharge of FL-BOND II

Durable bonding strengths

The incorporation of an enhanced adhesive monomer and new photo initiator gives FL-BOND II exceptionally balanced bonding to both enamel and dentin. The HEMA-free Primer prevents the deterioration of the bonding interface due to hydrolysis that ensures long-term durability of the bond.



¹ Gordan V.V., Mondragon E., Watson R.E., Garvan C., Mjör I. A.: A clinical evaluation of a selfetching primer and giomer restorative material. Result at eight years. J Am Dent Assoc 138: 621-627 (2007)

Order information



FL-BOND II Kit

- Primer 5 ml
- Bonding Agent 5 ml
- MicroBrush pink and yellow, 25 pcs each
- 1 Petit Dish, 1 Plastic Cover for Petit Dish

PN 1306

FL-BOND II

- Primer 5 ml

PN 1307

FL-BOND II

- Bonding Agent 5 ml

PN 1308

FL-BOND II

- Etchant 6 ml

PN 1309



The filling composite Beautifil II is coordinated with FL-BOND II and convinces by an esthetic appearance and excellent handling properties. Combined with FL-BOND II the fluoride release is reinforced.



BEAUTIFIL II 6 Color Set

PN 1775

BEAUTIFIL II Tips 6 Color Set

PN 1776

BEAUTIFIL II Cosmetic Kit

PN 1777

BEAUTIFIL II Tips 3 Color Set

PN 1778



SHOFU INC. 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983, Japan
SHOFU INC. Singapore Branch 10 Science Park Road #03-12, The Alpha, Science Park II, Singapore 117684
SHOFU DENTAL CORPORATION 1225 Stone Drive, San Marcos, CA 92078-4059, USA
SHOFU DENTAL PRODUCTS LTD. Duke's Factory, Chiddingstone Causeway, Tonbridge, Kent, TN11 8 JU, UK
SHOFU DENTAL SUPPLIES (Shanghai) CO. LTD. No.645 Jiye Road, Sheshan Industry Park, Songjiang 201602, Shanghai, China
SHOFU DENTAL GmbH Am Brüll 17, 40878 Ratingen, Germany, Tel.: +49 2102 8664-0, Fax: +49 2102 8664-65, E-Mail: info@shofu.de