

from 0 to 80 in 60 sec.

# the new clear class



## glassbite

**the first clear transparent bite registration material with 80 Shore A**

- **thixotropic + stable**, does not flow into interdental spaces
- **high final hardness**, no springiness, ideal for further working
- **extremely short time in the mouth**, sets fully in only 60 sec.
- with **contouring tips**: for shallow, economical application and precise recording of relevant work area

new

### versatile uses

- ... as **matrix material** during light curing of composites, positioning aid
- ... as stamp material, for light-curing fillings (inlays & onlays)
- ... for labial impression with light-curing veneers and shells
- ... in orthodontics, for fixing and positioning brackets
- ... as index in the veneering technique and when making temporaries
- ... as wax-up key

