

*Super fast set, maximum hardness & fresh apple flavour*

# the cutting edge



## greenbite *apple*

**bite registration material, A-silicone based**

- Super fast set: extremely short time in the mouth, **only 60 seconds.**
- **Maximum final hardness**, approx. 45 Shore D, allows easy processing (cut and shape).
- **12 contouring tips free** in each standard packing for precise and easy application directly onto the teeth. Thixotropic, does not flow into interdental spaces.
- With **apple flavour** for stress-free bite taking

THE DENTAL  
ADVISOR  
++++½



**EXTRA:**  
12 contouring tips

