

for perfectionists

*easy,
just better*



Implant Mask

flexible gingival mask - easy to trim


The morphology of the gingival is a precious source of information for designing high-class implants. Implant Mask has particularly been developed **for this application**.



Equal flow around the abutments due to excellent flow properties. Detail reproduction due to innovative filler matrix on the base of A-silicone.



Perfect reproduction: Gingival sulcus, interdental papilla and all important tissue sections. Unlimited re-application, permanent dimensional stability.



Short setting time. Perfect working (trimming, grinding) through optimum final hardness of 70 shore A.

**mini-mix
cartridge:**
free-hand
application

