

# an intelligent combination

*the safety grip  
impression*



## monoprint *supra*

transfer silicone - monophase

Taking an implant impression makes special demands of the impression material. monoprint supra was developed **specifically for this**. An intelligent combination of relevant material characteristics. Nothing else.

- Stable tray consistency with above-average fluidity, combined with high wettability.
- Secure post fixation through **high final hardness**, snap set, combined with excellent elasticity and the greatest dimensional accuracy, neutral flavour.
- Circular embedding of the impression posts through uniform flow, even in difficult situations.

**plus**  
5 ml  
adhesive-A

