

		TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT							
Diameter/ Platform/ Connection angle	Name of the implant	Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker								
20°	20° UniAbutment 3.5/4.0	BSAD0731		ZBAD0731	24893 (0.5 mm)	24894 (1 mm)	24895 (2 mm)	24896 (4 mm)	24897 (6 mm)	24389 (8 mm)		
	20° UniAbutment 4.5/5.0				24398 (0.5 mm)	24399 (1 mm)	24400 (2 mm)	24401 (4 mm)	24402 (6 mm)	24403 (8 mm)		
45°	45° UniAbutment 3.5/4.0	BSAD0732		ZBAD0732	24899 (0.5 mm)	24900 (1 mm)	24901 (2 mm)	24902 (4 mm)	24396 (6 mm)	24397 (8 mm)		
	45° UniAbutment 4.5/5.0				24406 (0.5 mm)	24407 (1 mm)	24408 (2 mm)	24409 (4 mm)	24410 (6 mm)	24411 (8 mm)		

Diameter/ Platform/ Connection angle	Name of the implant	TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT							
		Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker								
Ø 3.0	OsseoSpeed™ TX 3.0 S	BSAD0131	BSAB0181	ZBAD0131				24982 (11 mm)	24983 (13 mm)	24984 (15 mm)		
Ø 3.5 - 4.0	OsseoSpeed™ TX 3.5 S	BSAD0132	BSAB0182	ZBAD0132		24930 (8 mm)	24931 (9 mm)	24932 (11 mm)	24933 (13 mm)	24934 (15 mm)	24935 (17 mm)	
	OsseoSpeed™ TX 4.0 S				24939 (6 mm)	24940 (8 mm)	24941 (9 mm)	24942 (11 mm)	24943 (13 mm)	24944 (15 mm)	24944 (17 mm)	
Ø 4.5 - 5.0	OsseoSpeed™ TX 4.5	BSAD0133	BSAB0183	ZBAD0133			24951 (9 mm)	24952 (11 mm)	24953 (13 mm)	24954 (15 mm)	24955 (17 mm)	
	OsseoSpeed™ TX 5.0						24961 (9 mm)	24962 (11 mm)	24963 (13 mm)	24964 (15 mm)	24965 (17 mm)	
	OsseoSpeed™ TX 5.0 S						24971 (9 mm)	24972 (11 mm)	24973 (13 mm)	24974 (15 mm)	24975 (17 mm)	

Diameter/ Platform/ Connection angle	Name of the implant	TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT							
		Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker								
Ø 3.4	Zahnwurzelförmig Nanotite™ Certain® 3.25	BSAD0831	BSAB0881	ZBAD0831	NINT3285 (8.5 mm)	NINT3210 (10 mm)	NINT3211 (11.5 mm)	NINT3213 (13 mm)	NINT3215 (15 mm)			
	Parallelwandig Nanotite™ Certain® 3.25				NIOSM385 (8.5 mm)	NIOSM310 (10 mm)	NIOSM311 (11.5 mm)	NIOSM313 (13 mm)	NIOSM315 (15 mm)	NIOSM318 (18 mm)		
Ø 4.1	Zahnwurzelförmig Nanotite™ Certain® 4.00	BSAD0832	BSAB0882	ZBAD0832	NINT485 (8.5 mm)	NINT410 (10 mm)	NINT411 (11.5 mm)	NINT413 (13 mm)	NINT415 (15 mm)			
	Parallelwandig Nanotite™ Certain® 4.00				NIOSS485 (8.5 mm)	NIOSS410 (10 mm)	NIOSS411 (11.5 mm)	NIOSS413 (13 mm)	NIOSS415 (15 mm)	NIOSS418 (18 mm)		
Ø 5.0	Zahnwurzelförmig Nanotite™ Certain® 5.00	BSAD0833	BSAB0883	ZBAD0833	NINT585 (8.5 mm)	NINT510 (10 mm)	NINT511 (11.5 mm)	NINT513 (13 mm)	NINT515 (15 mm)			
	Parallelwandig Nanotite™ Certain® 5.00				NIOSS585 (8.5 mm)	NIOSS510 (10 mm)	NIOSS511 (11.5 mm)	NIOSS513 (13 mm)	NIOSS515 (15 mm)			
Ø 6.0	Zahnwurzelförmig Nanotite™ Certain® 6.00	BSAD0834	BSAB0884	ZBAD0834	NINT685 (8.5 mm)	NINT610 (10 mm)	NINT611 (11.5 mm)	NINT613 (13 mm)	NINT615 (15 mm)			
	Parallelwandig Nanotite™ Certain® 5.00				NIOSS685 (8.5 mm)	NIOSS610 (10 mm)	NIOSS611 (11.5 mm)	NIOSS613 (13 mm)	NIOSS615 (15 mm)			

Diameter/ Platform/ Connection angle	Name of the implant	TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT							
		Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker								
Ø 3.4	Parallelwandig Osseotite® 3.25	BSAD0931	BSAB0981	ZBAD0931		OSM385 (8.5 mm)	OSM310 (10 mm)	OSM311 (11 mm)	OSM313 (13 mm)	OSM315 (15 mm)	OSM318 (18 mm)	
Ø 4.1	Parallelwandig Osseotite® 3.75	BSAD0932	BSAB0982	ZBAD0932		OSS385 (8.5 mm)	OSS310 (10 mm)	OSS311 (11 mm)	OSS313 (13 mm)	OSS315 (15 mm)	OSS318 (18 mm)	OSS320 (20 mm)
	Parallelwandig Osseotite® 4.00					OSS485 (8.5 mm)	OSS410 (10 mm)	OSS411 (11 mm)	OSS413 (13 mm)	OSS415 (15 mm)	OSS418 (18 mm)	OSS420 (20 mm)
Ø 5.0	Parallelwandig Osseotite® 3.25	BSAD0933	BSAB0983	ZBAD0933	OSS507 (7 mm)	OSS585 (8.5 mm)	OSS510 (10 mm)	OSS511 (11 mm)	OSS513 (13 mm)	OSS515 (15 mm)	OSS518 (18 mm)	
Ø 6.0	Parallelwandig Osseotite® 3.25	BSAD0934	BSAB0984	ZBAD0934	OSS607 (7 mm)	OSS685 (8.5 mm)	OSS610 (10 mm)	OSS611 (11 mm)	OSS613 (13 mm)	OSS615 (15 mm)	OSS618 (18 mm)	

Diameter/ Platform/ Connection angle	Name of the implant	TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT							
		Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker								
Ø 3.3	ROOT-LINE Promote® 3.30	BSAD2031	BSAB2081	ZBAD2031								
	SCREW-LINE Promote® 3.30				K1042.3311 (11 mm)	K1042.3313 (13 mm)	K1042.3316 (16 mm)					
	ROOT-LINE Promote® plus 3.30				K1052.3311 (11 mm)	K1052.3313 (13 mm)	K1052.3316 (16 mm)					
Ø 3.8	ROOT-LINE Promote® 3.80	BSAD2032	BSAB2082	ZBAD2032	J1032.3809 (9 mm)	J1032.3811 (11 mm)	J1032.3813 (13 mm)	J1032.3816 (16 mm)				
	SCREW-LINE Promote® 3.80				K1042.3809 (9 mm)	K1042.3811 (11 mm)	K1042.3813 (13 mm)	K1042.3816 (16 mm)				
	ROOT-LINE Promote® plus 3.80				K1052.3809 (9 mm)	K1052.3811 (11 mm)	K1052.3813 (13 mm)	K1052.3816 (16 mm)				
	SCREW-LINE Guide Promote® plus 3.80				K1053.3809 (9 mm)	K1053.3811 (11 mm)	K1053.3813 (13 mm)					
Ø 4.3	ROOT-LINE Promote® 4.30	BSAD2033	BSAB2083	ZBAD2033	J1032.4309 (9 mm)	J1032.4311 (11 mm)	J1032.4313 (13 mm)	J1032.4316 (16 mm)				
	SCREW-LINE Promote® 4.30				K1042.4309 (9 mm)	K1042.4311 (11 mm)	K1042.4313 (13 mm)	K1042.4316 (16 mm)				
	ROOT-LINE Promote® plus 4.30				K1052.4309 (9 mm)	K1052.4311 (11 mm)	K1052.4313 (13 mm)	K1052.4316 (16 mm)				
	SCREW-LINE Guide Promote® plus 4.30				K1053.4309 (9 mm)	K1053.4311 (11 mm)	K1053.4313 (13 mm)					

Diameter/ Platform/ Connection angle	Name of the implant	TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT							
		Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker								
Ø 5.0	ROOT-LINE Promote® 5.00	BSAD2034	BSAB2084	ZBAD2034	J1032.5009 (9 mm)	J1032.5011 (11 mm)	J1032.5013 (13 mm)	J1032.5016 (16 mm)				
	SCREW-LINE Promote® 5.00				K1042.5009 (9 mm)	K1042.5011 (11 mm)	K1042.501 (13 mm)	K1042.5016 (16 mm)				
	ROOT-LINE Promote® plus 5.00				K1052.5009 (9 mm)	K1052.5011 (11 mm)	K1052.5013 (13 mm)	K1052.5016 (16 mm)				
Ø 6.0	ROOT-LINE Promote® 6.00	BSAD2035	BSAB2085	ZBAD2035	J1032.6009 (9 mm)	J1032.6011 (11 mm)	J1032.6013 (13 mm)	J1032.6016 (16 mm)				
	SCREW-LINE Promote® 6.00				K1042.6009 (9 mm)	K1042.6011 (11 mm)	K1042.6013 (13 mm)	K1042.6016 (16 mm)				
	ROOT-LINE Promote® plus 6.00				K1052.6009 (9 mm)	K1052.6011 (11 mm)	K1052.6013 (13 mm)	K1052.6016 (16 mm)				

Diameter/ Platform/ Connection angle	Name of the implant	TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT							
		Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker								
Ø 3.3 3.8 4.3	Stegaufbau 3.30	BSAD2131		ZBAD2131	J2255.3305 (0.5 mm)	J2255.3320 (2 mm)						
	Stegaufbau 3.80				J2255.3805 (0.5 mm)	J2255.3820 (2 mm)	J2255.3840 (4 mm)					
	Stegaufbau 4.30				J2255.4305 (0.5 mm)	J2255.4320 (2 mm)	J2255.4340 (4 mm)					
Ø 5.0 6.0	Stegaufbau 5.00	BSAD2132		ZBAD2132	J2255.5005 (0.5 mm)	J2255.5020 (2 mm)	J2255.5040 (4 mm)					
	Stegaufbau 5.00				J2255.6005 (0.5 mm)	J2255.6020 (2 mm)	J2255.6040 (4 mm)					

Friident DENTSPLY

XiVE®

Diameter/ Platform/ Connection angle	Name of the implant	TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT							
		Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker								
Ø 3.4	XiVE® S plus 3.40	BSAD1031	BSAB1081	ZBAD1031	26 – 2431 (9.5 mm)	26 – 2432 (11 mm)	26 – 2433 (13 mm)	26 – 2435 (15 mm)	26 – 2438 (18 mm)			
Ø 3.8	XiVE® S plus 3.80	BSAD1032	BSAB1082	ZBAD1032	26 – 2440 (8 mm)	26 – 2441 (9.5 mm)	26 – 2442 (11 mm)	26 – 2443 (13 mm)	26 – 2445 (15 mm)	26 – 2448 (18 mm)		
Ø 4.5	XiVE® S plus 4.5	BSAD1033	BSAB1083	ZBAD1033	26 – 2450 (8 mm)	26 – 2451 (9.5 mm)	26 – 2452 (11 mm)	26 – 2453 (13 mm)	26 – 2455 (15 mm)	26 – 2458 (18 mm)		
Ø 5.5	XiVE® S plus 5.5	BSAD1034	BSAB1084	ZBAD1034	26 – 2460 (8 mm)	26 – 2461 (9.5 mm)	26 – 2462 (11 mm)	26 – 2463 (13 mm)	26 – 2465 (15 mm)			

		TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT								
Diameter/ Platform/ Connection angle	Name of the implant	Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker									
Ø 3.5	MK III Groovy NP (3.30)	BSAD0231	BSAB0281	ZBAD0231		32114 (10 mm)	32115 (11.5 mm)	32116 (13 mm)	32117 (15 mm)				
	MK III TiUnite® NP (3.30)				28906 (10 mm)	28907 (11.5 mm)	28908 (13 mm)	28909 (15 mm)					
	NobelSpeedy™ Groovy NP (3.30)				33123 (10 mm)	33124 (11.5 mm)	33125 (13 mm)	33126 (15 mm)					
Ø 4.1	Mk III Shorty RP (3.75)	BSAD0232	BSAB0282	ZBAD0232	32118 (7 mm)								
	MK III Groovy RP (3.75)					32119 (8.5 mm)	32120 (10 mm)	32121 (11.5 mm)	32122 (13 mm)	32123 (15 mm)	32124 (18 mm)		
	Mk III TiUnite® RP (3.75)				28911 (7 mm)	28912 (8.5 mm)	28913 (10 mm)	28914 (11.5 mm)	28915 (13 mm)	28916 (15 mm)	28917 (18 mm)		
	Mk III Shorty RP (4.00)				32125 (7 mm)								
	MK III Groovy RP (4.00)					32126 (8.5 mm)	32127 (10 mm)	32128 (11.5 mm)	32129 (13 mm)	32130 (15 mm)	32131 (18 mm)		
	Mk III TiUnite® RP (4.00)				28918 (7 mm)	28919 (8.5 mm)	28920 (10 mm)	28921 (11.5 mm)	28922 (13 mm)	28923 (15 mm)	28924 (18 mm)		
	NobelSpeedy™ Shorty RP (4.00)				32146 (7 mm)								
	NobelSpeedy™ Groovy RP (4.00)					32147 (8.5 mm)	32148 (10 mm)	32149 (11.5 mm)	32150 (13 mm)	32151 (15 mm)	32152 (18 mm)		
MK IV TiUnite® RP (4.00)	28932 (7 mm)	28933 (8.5 mm)	28934 (10 mm)	28935 (11.5 mm)	28936 (13 mm)	28937 (15 mm)	28938 (18 mm)						

		TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT							
Diameter/ Platform/ Connection angle	Name of the implant	Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker								
Ø 5.1	Mk III Shorty WP (5.00)	BSAD0233	BSAB0283	ZBAD0233	32132 (7 mm)							
	Mk III Groovy WP (5.00)					32133 (8.5 mm)	32134 (10 mm)	32135 (11.5 mm)	32136 (13 mm)	32137 (15 mm)	32138 (18 mm)	
	Mk III TiUnite® WP (5.00)				28925 (7 mm)	28926 (8.5 mm)	28927 (10 mm)	28928 (11.5 mm)	28929 (13 mm)	28930 (15 mm)		
	NobelSpeedy™ Shorty WP (5.00)				32153 (7 mm)							
	NobelSpeedy™ Groovy WP (5.00)					32154 (8.5 mm)	32155 (10 mm)	32156 (11.5 mm)	32157 (13 mm)	32158 (15 mm)	32159 (18 mm)	
	Mk IV TiUnite® WP (5.00)				28939 (7 mm)	28940 (8.5 mm)	28941 (10 mm)	28942 (11.5 mm)	28943 (13 mm)	28944 (15 mm)		
	NobelSpeedy™ Shorty 6.00				32139 (7 mm)							
	NobelSpeedy™ Groovy 6.00					32140 (8.5 mm)	32141 (10 mm)	32142 (11.5 mm)	32143 (13 mm)	32144 (15 mm)	32145 (18 mm)	

Diameter/ Platform/ Connection angle	Name of the implant	TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT								
		Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker									
NP-RP	NP Brånemark	BSAD0431		ZBAD0431	29176 (1 mm)	29177 (2 mm)	29178 (3 mm)						
	RP Brånemark				29179 (1 mm)	29180 (2 mm)	29181 (3 mm)	29182 (4 mm)	29184 (5 mm)				
	NP Active				34186 (1.5 mm)	34187 (2.5 mm)	34595 (3.5 mm)						
	RP Active				24190 (1.5 mm)	34191 (2.5 mm)	34596 (3.5 mm)	35407 (4.5 mm)					
	NP Replace				29196 (1 mm)	29197 (2 mm)	29198 (3 mm)						
	RP Replace				29199 (1 mm)	29200 (2 mm)	29201 (3 mm)	29202 (4 mm)	29203 (5 mm)				
WP	WP Brånemark	BSAD0432		ZBAD0432	29184 (1 mm)	29185 (2 mm)	29186 (3 mm)						
	WP Replace				29204 (1 mm)	29205 (2 mm)	29206 (3 mm)						

		TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT							
Diameter/ Platform/ Connection angle	Name of the implant	Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker								
Ø 3.5	Tapered NP (3.50)	BSAD0331	BSAB0381	ZBAD0331	32211 (8 mm)	32212 (10 mm)	36100 (11.5 mm)	32213 (13 mm)	32214 (16 mm)			
	Select Tapered NP (3.50)				36104 (8 mm)	29401 (10 mm)	36105 (11.5 mm)	29402 (13 mm)	29403 (16 mm)			
	Straight NP (3.50)					32160 (10 mm)	32161 (11.5 mm)	32162 (13 mm)	32163 (15 mm)			
	Select Straight NP (3.50)					28946 (10 mm)	28947 (11.5 mm)	28948 (13 mm)	28949 (15 mm)			
	NobelSpeedy Replace® NP (3.50)					32196 (10 mm)	32197 (11.5 mm)	32198 (13 mm)	32199 (15 mm)			
	Select TC NP (3.50)					36285 (7+3 mm)	36286 (10+3 mm)	36287 (13+3 mm)	36288 (15+3 mm)			
Ø 4.3	Tapered RP (4.30)	BSAD0332	BSAB0382	ZBAD0332	32215 (8 mm)	32216 (10 mm)	36101 (11.5 mm)	32217 (13 mm)	32218 (16 mm)			
	Select Tapered RP (4.30)				36106 (8 mm)	29413 (10 mm)	36107 (11.5 mm)	29414 (13 mm)	29415 (16 mm)			
	Straight RP (4.30)					32181 (10 mm)	32182 (11.5 mm)	32183 (13 mm)	32184 (15 mm)	32185 (18 mm)		
	Select Straight RP (4.30)					28951 (10 mm)	28952 (11.5 mm)	28953 (13 mm)	28954 (15 mm)	28955 (18 mm)		
	NobelSpeedy Replace RP (4.30)					32200 (10 mm)	32201 (11.5 mm)	32202 (13 mm)	32203 (15 mm)	32204 (18 mm)		
	Select TC RP (4.30)					36289 (7+3 mm)	36290 (10+3 mm)	36291 (13+3 mm)	36292 (15+3 mm)			

Nobel Biocare®

NobelReplace®

		TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT							
Diameter/ Platform/ Connection angle	Name of the implant	Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker								
Ø 5.0	Tapered WP (5.00)	BSAD0333	BSAB0383	ZBAD0333	32219 (8 mm)	32220 (10 mm)	36102 (11.5 mm)	32221 (13 mm)	32222 (16 mm)			
	Select Tapered WP (5.00)				36108 (8 mm)	29423 (10 mm)	36109 (11.5 mm)	29424 (13 mm)	29425 (16 mm)			
	Straight WP (5.00)					32186 (10 mm)	32187 (11.5 mm)	32188 (13 mm)	32189 (15 mm)	32190 (18 mm)		
	Select Straight WP (5.00)					28956 (10 mm)	28957 (11.5 mm)	28958 (13 mm)	28959 (15 mm)			
	NobelSpeedy Replace WP (5.00)					32205 (10 mm)	32206 (11.5 mm)	32207 (13 mm)	32208 (15 mm)	32209 (18 mm)		
Ø 6.0	Tapered 6.00	BSAD0334	BSAB0384	ZBAD0334	32223 (8 mm)	32224 (10 mm)	36103 (11.5 mm)	32225 (13 mm)	32226 (16 mm)			
	Select Tapered 6.00				36110 (8 mm)	32949 (10 mm)	36111 (11.5 mm)	32950 (13 mm)				
	NobelSpeedy Replace 6.00					32191 (10 mm)	32192 (11.5 mm)	32193 (13 mm)	32194 (15 mm)	32195 (18 mm)		

		TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT							
Diameter/ Platform/ Connection angle	Name of the implant	Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker								
Ø 3.3	Sicace 3.40	BSAD1231	BSAB1281	ZBAD1231	935182 (7.5 mm)	935170 (9.5 mm)	935171 (11.5 mm)	935172 (13 mm)	935173 (14.5 mm)			
	Sicace 4.00				935183 (7.5 mm)	935174 (9.5 mm)	935175 (11.5 mm)	935176 (13 mm)	935177 (14.5 mm)			
Ø 4.2	Sicace 5.00	BSAD1232	BSAB1282	ZBAD1232	935184 (7.5 mm)	935178 (9.5 mm)	935179 (11.5 mm)	935180 (13 mm)	935181 (14.5 mm)			

Diameter/ Platform/ Connection angle	Name of the implant	TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT							
		Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker								
Ø 3.3	SLActive BL 3.3 NC	BSAD0531	BSAB0581	ZBAD0531	021.2108 (8 mm)	021.2110 (10 mm)	021.2112 (12 mm)	021.2114 (14 mm)				
Ø 4.1	SLActive BL 4.1 RC	BSAD0532	BSAB0582	ZBAD0532	021.4108 (8 mm)	021.4110 (10 mm)	021.4112 (12 mm)	021.4114 (14 mm)				
	SLActive BL 4.8 RC				021.6108 (8 mm)	021.6110 (10 mm)	021.6112 (12 mm)	021.6114 (14 mm)				

		TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT								
Diameter/ Platform/ Connection angle	Name of the implant	Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker									
NN Ø 3.5	SLActive SP 3.3 NN	BSAD0631	BSAB0681	ZBAD0631		033.951S (8 mm)	033.952S (10 mm)	033.953S (12 mm)	033.954S (14 mm)				
RN Ø 4.8	SLActive SP 3.3 RN	BSAD0632	BSAB0682	ZBAD0632 (direkt auf Implantat) ZBAD0642 (mit synOcta-Adapter)		033.151S (8 mm)	033.152S (10 mm)	033.153S (12 mm)	033.154S (14 mm)				
	SLActive SP 4.1 RN				033.050S (6 mm)	033.051S (8 mm)	033.052S (10 mm)	033.053S (12 mm)	033.054S (14 mm)				
	SLActive SP 4.8 RN				033.250S (6 mm)	033.251S (8 mm)	033.252S (10 mm)	033.253S (12 mm)	033.254S (14 mm)				
WN Ø 6.5	SLActive SP 4.8 WN	BSAD0633	BSAB0683	ZBAD0633 (direkt auf Implantat) ZBAD0643 (mit synOcta-Adapter)	033.650S (6 mm)	033.651S (8 mm)	033.652S (10 mm)	033.653S (12 mm)					

Sweden & Martina

Kohno/Premium

		TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT							
Diameter/ Platform/ Connection angle	Name of the implant	Set CONICAL cemented titanium base NON HEX	Set PARALLELL cemented titanium base HEX	Set Scanmarker								
Ø 3.5	Premium Straight ZirTi 3.30	BSAD1531		ZBAD1531		A-ZT-330-100 (10 mm)	A-ZT-330-115 (11.5 mm)	A-ZT-330-130 (13 mm)	A-ZT-330-150 (15 mm)			
	Premium Straight DES 3.30				A-DS-330-100 (10 mm)	A-DS-330-115 (11.5 mm)	A-DS-330-130 (13 mm)	A-DS-330-150 (15 mm)				
Ø 3.8	Premium SP ZirTi 3.80	BSAD1532		ZBAD1532	A-ZT-380SP-085 (8.5 mm)	A-ZT-380SP-100 (10 mm)	A-ZT-380SP-115 (11.5 mm)	A-ZT-380SP-130 (13 mm)	A-ZT-380SP-150 (15 mm)			
	Premium SP DES 3.80		A-DS-380SP-085 (8.5 mm)		A-DS-380SP-100 (10 mm)	A-DS-380SP-115 (11.5 mm)	A-DS-380SP-130 (13 mm)	A-DS-380SP-150 (15 mm)				
	Premium Straight ZirTi 3.80		A-ZT-380-085 (8.5 mm)		A-ZT-380-100 (10 mm)	A-ZT-380-115 (11.5 mm)	A-ZT-380-130 (13 mm)	A-ZT-380-150 (15 mm)	A-ZT-380-180 (18 mm)			
	Premium Straight DES 3.80		A-DS-380-085 (8.5 mm)		A-DS-380-100 (10 mm)	A-DS-380-115 (11.5 mm)	A-DS-380-130 (13 mm)	A-DS-380-150 (15 mm)				
	Kohno SP DES 3.80				K-DS-380SP-100 (10 mm)	K-DS-380SP-115 (11.5 mm)	K-DS-380SP-130 (13 mm)	K-DS-380SP-150 (15 mm)				
	Kohno Straight DES 3.80				K-ZT-380-100 (10 mm)	K-ZT-380-115 (11.5 mm)	K-ZT-380-130 (13 mm)	K-ZT-380-150 (15 mm)	K-ZT-380-180 (18 mm)			
Ø 4.25	Premium SP ZirTi 4.25	BSAD1533		ZBAD1533	A-ZT-425SP-085 (8.5 mm)	A-ZT-425SP-100 (10 mm)	A-ZT-425SP-115 (11.5 mm)	A-ZT-425SP-130 (13 mm)	A-ZT-425SP-150 (15 mm)			
	Premium SP DES 4.25		A-DS-425SP-085 (8.5 mm)		A-DS-425SP-100 (10 mm)	A-DS-425SP-115 (11.5 mm)	A-DS-425SP-130 (13 mm)	A-DS-425SP-150 (15 mm)				
	Premium Straight ZirTi 4.25		A-ZT-425-085 (8.5 mm)		A-ZT-425-100 (10 mm)	A-ZT-425-115 (11.5 mm)	A-ZT-425-130 (13 mm)	A-ZT-425-150 (15 mm)	A-ZT-425-180 (18 mm)			
	Premium Straight DES 4.25		A-DS-425-085 (8.5 mm)		A-DS-425-100 (10 mm)	A-DS-425-115 (11.5 mm)	A-DS-425-130 (13 mm)	A-DS-425-150 (15 mm)				
	Kohno SP DES 4.25				K-DS-425SP-100 (10 mm)	K-DS-425SP-115 (11.5 mm)	K-DS-425SP-130 (13 mm)	K-DS-425SP-150 (15 mm)				
	Kohno Straight DES 4.25				K-DS-425-100 (10 mm)	K-DS-425-115 (11.5 mm)	K-DS-425-130 (13 mm)	K-DS-425-150 (15 mm)	K-DS-425-180 (18 mm)			

Sweden & Martina

Kohno/Premium

		TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT								
Diameter/ Platform/ Connection angle	Name of the implant	Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker									
Ø 5.0	Premium SP ZirTi 5.00	BSAD1534		ZBAD1534	A-ZT-500SP-085 (8.5 mm)	A-ZT-500SP-100 (10 mm)	A-ZT-500SP-115 (11.5 mm)	A-ZT-500SP-130 (13 mm)	A-ZT-500SP-150 (15 mm)				
	Premium SP DES 5.00				A-DS-500SP-085 (8.5 mm)	A-DS-500SP-100 (10 mm)	A-DS-500SP-115 (11.5 mm)	A-DS-500SP-130 (13 mm)	A-DS-500SP-150 (15 mm)				
	Premium Straight ZirTi 5.00				A-ZT-500-085 (8.5 mm)	A-ZT-500-100 (10 mm)	A-ZT-500-115 (11.5 mm)	A-ZT-500-130 (13 mm)	A-ZT-500-150 (15 mm)				
	Premium Straight DES 5.00				A-DS-500-085 (8.5 mm)	A-DS-500-100 (10 mm)	A-DS-500-115 (11.5 mm)	A-DS-500-130 (13 mm)	A-DS-500-150 (15 mm)				
	Kohno SP DES 5.00					K-DS-500SP-100 (10 mm)	K-DS-500SP-115 (11.5 mm)	K-DS-500SP-130 (13 mm)					
	Kohno Straight DES 5.00					K-DS-500-100 (10 mm)	K-DS-500-115 (11.5 mm)	K-DS-500-130 (13 mm)					
Ø 6.0		BSAD1535		ZBAD1535	K-DS-600-085 (8.5 mm)	K-DS-600-100 (10 mm)	K-DS-600-115 (11.5 mm)						

Sweden & Martina

Multi Unit Abutment

		TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT							
Diameter/ Platform/ Connection angle	Name of the implant	Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker								
MUA	Plattform 3.30	BSAD1431		ZBAD1431	A-PAD-AD330-15	A-PAD-AD330-30	A-PAD-AD330-40					
	Plattform 3.80				A-PAD-AD380-15	A-PAD-AD380-30	A-PAD-AD380-40					
	Plattform 4.25				A-PAD-AD425-15	A-PAD-AD425-30	A-PAD-AD425-40					
	Plattform 5.00				A-PAD-AD500-15	A-PAD-AD500-30	A-PAD-AD500-40					

Sweden & Martina

Out-Link

		TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT								
Diameter/ Platform/ Connection angle	Name of the implant	Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker									
Ø 3.3	Outlink 3S 3.30	BSAD1331		ZBAD1331				E-3S-330-100 (10 mm)	E-3S-330-115 (11.5 mm)	E-3S-330-130 (13 mm)	E-3S-330-150 (15 mm)		
	Outlink ZirTi 3.30								E-3S-330-100 (10 mm)	E-3S-330-115 (11.5 mm)	E-3S-330-130 (13 mm)	E-3S-330-150 (15 mm)	
Ø 4.1	Outlink 3S 3.75	BSAD1332		ZBAD1332			E-3S-375-085 (8.5 mm)	E-3S-375-100 (10 mm)	E-3S-375-115 (11.5 mm)	E-3S-375-130 (13 mm)	E-3S-375-150 (15 mm)		
	Outlink ZirTi 3.75								E-ZT-375-085 (8.5 mm)	E-ZT-375-100 (10 mm)	E-ZT-375-115 (11.5 mm)	E-ZT-375-130 (13 mm)	E-ZT-375-150 (15 mm)
	Outlink 3S 4.10								E-3S-410-085 (8.5 mm)	E-3S-410-100 (10 mm)	E-3S-410-115 (11.5 mm)	E-3S-410-130 (13 mm)	E-3S-410-150 (15 mm)
	Outlink ZirTi 4.10								E-ZT-410-085 (8.5 mm)	E-ZT-410-100 (10 mm)	E-ZT-410-115 (11.5 mm)	E-ZT-410-130 (13 mm)	E-ZT-410-150 (15 mm)
Ø 5.0	Outlink 3S 5.00	BSAD1333		ZBAD1333			E-3S-500-085 (8.5 mm)	E-3S-500-100 (10 mm)	E-3S-500-115 (11.5 mm)	E-3S-500-130 (13 mm)			
	Outlink ZirTi 5.00								E-ZT-500-050 (5 mm)	E-ZT-500-070 (7 mm)	E-ZT-500-085 (8.5 mm)	E-ZT-500-100 (10 mm)	E-ZT-500-115 (11.5 mm)

		TITANIUM BASES AND ACCESSORIES – ITEM NUMBER ZIRKONZAHN			ITEM NUMBER ORIGINAL IMPLANT							
Diameter/ Platform/ Connection angle	Name of the implant	Set CONICAL cemented titanium base NON HEX	Set PARALELL cemented titanium base HEX	Set Scanmarker								
Ø 3.5	MTX 3.70	BSAD1131	BSAB1181	ZBAD1131	TSVB8 (8 mm)	TSVB10 (10 mm)	TSVB11 (11.5 mm)	TSVB13 (13 mm)	TSVB16 (16 mm)			
	MP-1*HA 3.70				TSVH8 (8 mm)	TSVH10 (10 mm)	TSVH11 (11.5 mm)	TSVH13 (13 mm)	TSVH16 (16 mm)			
	MTX 4.10				TSV4B8 (8 mm)	TSV4B10 (10 mm)	TSV4B11 (11.5 mm)	TSV4B13 (13 mm)	TSV4B16 (16 mm)			
	MP-1*HA 4.10				TSV4H8 (8 mm)	TSV4H10 (10 mm)	TSV4H11 (11.5 mm)	TSV4H13 (13 mm)	TSV4H16 (16 mm)			
Ø 4.5	MTX 4.70	BSAD1132	BSAB1182	ZBAD1132	TSVWB8 (8 mm)	TSVWB10 (10 mm)	TSVWB11 (11.5 mm)	TSVWB13 (13 mm)	TSVWB16 (16 mm)			
	MP-1*HA 4.70				TSVWH8 (8 mm)	TSVWH10 (10 mm)	TSVWH11 (11.5 mm)	TSVWH13 (13 mm)	TSVWH16 (16 mm)			
Ø 5.7	MTX 6.00	BSAD1133	BSAB1183	ZBAD1133	TSV6B8 (8 mm)	TSV6B10 (10 mm)	TSV6B11 (11.5 mm)	TSV6B13 (13 mm)	TSV6B16 (16 mm)			
	MP-1*HA 6.00				TSV6H8 (8 mm)	TSV6H10 (10 mm)	TSV6H11 (11.5 mm)	TSV6H13 (13 mm)	TSV6H16 (16 mm)			