


NEW


MILLING UNIT M1 ABUTMENT

An overview of the special features

- For milling of titan bases and abutments made of pre-fabricated titanium abutment blanks
- 4+1 axes simultaneous milling technology
- Wet processing
- Milling Spindle Hard Automatic with fivefold tool changer
- Space-saving
- Modular
- Equipped with sealing air
- Processable materials: pre-milled abutment blanks made of titanium, glass ceramics
- 47.7 x 69.3 x 61.3 cm
- Weight: about 80 kg



4 AXES



MODULAR



TOOLCHANGER
5x



WET



METAL



PRE-MILL
1x

All information is subject to change. Errors and omissions excepted. Version: 07.03.2013

