



Smile Line ☺  
made in switzerland

## AGATE DRY MIXING TRAY

THE AUTHENTIC BEAUTY OF THE NATURE

Your pre-mixed materials (paste-opaques, pre-mixed stains) are always ready to be used thanks to the new Smile Line Agate dry mixing trays.

The Agate, a semi-precious stone,

is cut, transformed and machined manually. Its uniqueness and beauty are enhanced in the ergonomic Smile Line translucent container.

Because of the nature of the material,

every single piece is a surprise for the user. Smile Line proposes one tray with 16 wells, intended for your paste-opaques and another with 20 wells, for your stains in paste.

### DESCRIPTION

AGATE dry mixing tray for stains in paste, 20 wells  
AGATE dry mixing tray for paste-opaques, 16 wells

### N° Art.

1045  
1055

IMPORTER:

mail@smileline.ch www.smileline.ch

Smile Line ☺  
made in switzerland