


**Zirkonzahn<sup>®</sup>**

*Human Zirconium Technology*

**NEW**

**HEROES COLLECTION**

*by Enrico Steger*





*Cut-Back FIRE*



AIDA



AREUS



DEMI



HERMES



PAIAN





# WONDERS OF NATURE

*Teeth, formed by nature, appear in the most beautiful shapes.*


*Among thousand models and with meticulous accuracy I've selected these 10 marvels, which I would now like to share with you.*

*We've made 10 books and we've implemented all teeth into the CAD/CAM software.*

*Luís Stejer*

## THE VIRTUAL AND INTELLIGENT TOOTH LIBRARY

- 10 natural and aesthetically pleasing tooth sets
- Fully anatomic or in 4 cut-back designs FIRE, WATER, AIR, EARTH for veneering with ceramics
- Aesthetic base for any kind of restoration: single crowns, small bridges, fully anatomic bridges (Prettau® Bridges), complete dentures
- Individual adjustments for any patient possible
- Combinable with the natural tooth substance, e.g. if only the occlusal surface of the tooth is abraded
- Different occlusal concepts can be realised, e.g. canine or posterior tooth guidance
- Suitable for treatment planning



RHEIA



SAO




SERAPHIM





THALIA



TRITON



# Zirkonzahn®



## HUMAN ZIRCONIUM TECHNOLOGY

ZIRKONZAHN Worldwide – Tel +39 0474 066 680 – info@zirkonzahn.com – www.zirkonzahn.com

THE HEROES COLLECTION  
IS ALSO AVAILABLE IN  
BOOK FORMAT!

FURHER INFORMATION:  
[WWW.ZIRKONZAHN.COM](http://WWW.ZIRKONZAHN.COM)



WEAA7721=

All information is subject to change. Errors and omissions excepted. Version: 27.09.2013