

IN-VITRO COMPARATIVE STUDY ON FLEXURAL STRENGTH OF PRETTAU® ANTERIOR®


University of Pennsylvania April 2016

In-vitro comparative study of Prettau® Anterior® Zirconia


ABSTRACT

Scientifically proven – 10-year stability of Prettau® Anterior® . No decrease of flexural strength after 10 years (2,4 million cycles).

An in-vitro comparative study performed by the University of Pennsylvania has shown that there is no decrease in the flexural strength of Prettau® Anterior® after 2,4 Millions cycles of thermo-mechanical chewing simulation, which is equivalent to 10 years of clinical use.



Markus B. Blatz,
DMD, PhD Philadelphia PA,



Michael M. Bergler, CBE, MDT
Philadelphia PA,