

FIBER KOČIĆI sa ResiCem SHOFU

Dr Darko Stamatović

IKODENTAL
DENTAL DEPO



ResiCem twinlock cement i Porcelan primer



Preparacije korena kanala



Kap ResiCem bonda primer A



Kap u isto udubljenje ResiCem bonda primer B i pomešati



Premazati unutrašnjost preparisanog kanala



U drugo udubljenje na pločici kapnuti kap Porcelan primera



Četkicom premazati fiber kočić i ostaviti 10 sek.



Iz tube sa automixerom direktno u kanal aplikovati twinlock cement 1



Iz tube sa automixerom direktno u kanal aplikovati twinlock cement 2



Odmah ubaciti fiber kočiće



Odmah ubaciti fiber kočiće, prosvetlite 20 sek. i sačekati 4min da se cement skroz veže



Posle 4 minuta



Izgled bosle brušenja



Ro snimak pre postavljanja post kočića



Ro snimak posle postavljanja post kočića

* Zahvaljujemo Dr Darku Stamatoviću na podršci

ResiCem®

Dual-stvrdnjavajući, univerzalni sistem za cementiranje, idealan za alumina i cirkonija keramiku.

Rendgen kontrastan, dual-stvrdnjavajući vezivni sistem sa samokondicionirajućim prajmerom A i B za adhezivno cementiranje indirektnih restauracija svih vrsta keramika, kompozita, plemenitih ili neplemenitih legura. Odlikuju ga izvanredne fizičke osobine, koje omogućavaju pouzdano cementiranje, jednostavnu upotrebu.

ResiCem je depo fluorida, jer sadrži PRG – Shofu patentirani punjač, koji obezbeđuje antikariogenu zaštitu. Proverite zašto je jedan od najtraženijih Shofu proizvoda ! Dostupan u bojama: providnoj(clear) i slonovača (Ivory), idealan za sve estetske zahteve, pri apsolutnoj stabilnosti boje.



Odmah nakon cementiranja



2 min. nakon cementiranja



PN/šifra: 3220

ResiCem Basic set: ResiCem pasta „Clear“, auto-mixing, 2 x 5 ml, Primer A 3 ml, Primer B 2.5 ml, Microtip fine (1 drška, 50 vrhova), 20 kanila za mešanje, posuda za mešanje.

IKODENTAL
DENTAL DEPO

